


RÈGLES ÉLASTOMÈRES AIMANTÉES

Arase

Longueur 1.25 ml



A (mm)	B (mm)	C (mm)
36	30	30
40	30	15
40	30	20




Architectonique

Longueur 1.25 ml



A (mm)	B (mm)	C (mm)
30	20	20
30	22	15
40	30	20



LITEAU ÉLASTOMÈRE AIMANTÉ

Aimant (1 côté)

Longueur 1.25 ml



A (mm)	B (mm)
15	15
20	20
30	30

